

PD Zinc Phosphate Cement C & B

Zinc Phosphate Cement

Composition

Powder: Zinc Oxide 88%, Magnesium Oxide 10%, Alu. Phosphate 2%

Liquid: Phosphoric Acid 58%, Alu. Hydroxide 7%, excip. ad 100%

Properties

An universal, all-purpose cement for cementing crowns, bridges, inlays, orthodontic appliances, etc. Also as filling material, a cavity liner & as a pulp capping under gold, amalgam, porcelain & acrylic fillings. Perfectly balanced, physically as well as chemically. It adheres just as strongly to natural dentines as to porcelain or metal & is resistant against buccal acids.

Packaging

30g powder, 15ml liquid

